

---

---

# Table of Contents

<b>Preface</b>	vii
<b>Foreword</b>	ix
<b>List of Figures</b>	xv
<b>List of Contributors</b>	xix
<b>A Principles of Biosafety</b>	1
A.1 Categories of Biosafety	2
A.2 Biohazard Materials	3
A.3 Recommended Work Practices	3
A.4 General Hygiene	4
A.4.1 Pipettes	5
A.4.2 Syringes and Scalpels	5
A.4.3 Biosafety Cabinets	5
A.4.4 Chemicals	6
A.4.5 Ultraviolet Light	7
A.4.6 General Housekeeping	7
A.5 Handling Biological Specimens	8
<b>B Introduction to the Cell Culture Laboratory</b>	11
B.1 The Laminar Flow Workbench — Biosafety Cabinet	12
B.2 The Centrifuge	12
B.3 The Waterbath	15
B.4 The Sink/Wet space	16
B.5 The Autoclave	16

B.6 The Refrigerator/Freezer	18
B.7 The Incubator	19
B.8 The Microscope	21
<b>C Protocols for Human Cell Culture</b>	<b>23</b>
C.1 Cell Isolation	25
C.2 Mesenchymal Cells	26
C.2.1 Adipose Cells	27
2.1.1 Enzymatic Digestion	27
C.2.2 Chondrocytes	28
2.2.1 Enzymatic Digestion	29
C.2.3 Endothelial Cells	30
2.3.1 HUVEC (human umbilical vein endothelial cells) Isolation	31
C.2.4 Fibroblasts	32
2.4.1 Explant Method	32
2.4.2 Enzymatic Digestion	33
2.4.3 Subculturing	34
C.2.5 Myocytes	34
2.5.1 Skeletal Muscle Cells	34
Enzymatic Digestion	35
Subculturing	35
2.5.2 Smooth Muscle Cells	36
Enzymatic Digestion	36
Subculturing	37
C.2.6 Osteoblasts	38
2.6.1 Wash-out Method	39
2.6.2 Explant Method	40
2.6.3 Enzymatic Digestion	41
2.6.4 Osteoblasts	41
2.6.5 Periosteal Cells	42
C.2.7 Dental Cells (Periodontal Ligament Fibroblasts & Alveolar Osteoblasts)	43
2.7.1 Explant Culture	43
Notes	46
C.3. Epithelial Cells	49
C.3.1 Respiratory Epithelial Cells	49
3.1.1 Enzymatic Digestion	50
C.3.2 Oral Mucosa Cells	51
3.2.1 Explant Method	51

3.2.2	Subculturing	52
C.3.3	Keratinocytes	54
3.3.1	Separating the Epidermis from Dermis	55
3.3.2	Isolating Single Cells from the Epidermis	55
3.3.3	Culturing on a Fibroblast Feeder Layer	55
3.3.4	Culturing in Serum-free Medium	56
3.3.5	Subculturing for both Feeder-layer and Serum-free	57
C.3.4	Hepatocytes	58
3.4.1	Enzymatic Digestion	58
C.3.5	Kidney Cells	60
3.5.1	Renal Cortical Cell Culture	60
3.5.2	Proximal Tubuli Cell Culture	61
3.5.3	Glomerula Mesangial Cell Culture	62
Notes		64
C.4	Bone Marrow Derived Cells	67
C.4.1	Bone Marrow Stromal Cells	67
C.4.2	Peripheral Blood Mononuclear Cells	69
Notes		72
<b>D</b>	<b>Reagents and Techniques for Cell Culture</b>	<b>75</b>
D.1	Selection of Cell Culture Media	76
D.1.1	Classical Media	76
D.1.2	Serum-Free and Speciality Media	77
Notes		83
D.2	Cell Culture Reagents and Media Supplements	86
D.2.1	Sera	86
2.1.1	Protocol for Isolation of Serum from a Blood Sample	86
D.2.2	Antibiotics (AB) and Antimycotics (AM)	87
D.2.3	Biological Buffers	88
D.2.4	Supplements	91
2.4.1	Water	91
2.4.2	L-Glutamine	91
2.4.3	Lipids	92
Notes		93
D.3	Cell Dissociation Products	96
D.3.1	General Considerations	96
D.3.2	Basic Strategy	97

D.3.3 Trypsin	98
D.3.4 Other Proteases	99
Notes	102
D.4 Quantitative Cell Proliferation Assays	105
D.4.1 MTS Cell Metabolic Assay	105
D.4.2 AlamarBlue™ Cell Metabolic Assay	106
D.4.3 WST-1™ Proliferation Assay	107
D.4.4 [ <sup>3</sup> H]-Thymidine Labelling Proliferation Assay	108
D.4.5 Pico-Green® DNA Quantification Assay	109
D.4.6 Use of Cell Proliferation Assays for High Cell Numbers	109
D.4.7 Colony Forming Unit Fibroblast Assay (CFU-F)	110
D.4.8 Growth Kinetics of Cell Populations	111
Notes	113
D.5 Transfection	116
D.5.1 Lipofection	116
5.1.1 Transfection with LIPOFECTAMINE™ (Invitrogen Life Technologies)	116
5.1.2 Transfection Optimisation	118
D.5.2 Electroporation	119
D.5.3 Calcium Phosphate Precipitation	119
D.5.4 Viral Vectors	119
Notes	121
D.6 Reference Guide	124
D.6.1 Troubleshooting	124
D.6.2 Cryopreservation and Thawing of Cells	126
D.6.3 Cell Counting with the Haemocytometer	130
Notes	133
Suggested Readings	137
Abbreviations	139
Milestones in the history of cell culture	141
Index	143